



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

StrepMan Magnet for 15 ml and 50 ml tubes

Cat. No.: 6-5650-065

Version: 2.4
Revision Date: 04.10.2021

Description	The StrepMan consists of a plastic frame with two embedded permanent magnets and accommodates two different tube sizes (for 2x 15 ml and 2x 50 ml conical reaction tubes). It can be used for magnetic cell isolation with microbeads such as Strep-Tactin® Magnetic Microbeads.
Attention	<p>Strong permanent magnet.</p> <p>The Strep-Tactin® Magnetic Microbeads (MB) have to be stored distant from strong magnetic fields like the StrepMan Magnet (at least 10 cm).</p> <p>Danger of injury by strong magnetic attraction of ferromagnetic materials such as iron.</p> <ul style="list-style-type: none">• Keep away from magnetic storage devices such as hard disks, diskettes, ID cards, credit cards, otherwise stored data may be deleted.• Keep away from electronic devices.• Persons with cardiac pacemakers should not be in close contact with this magnet.
Cleaning	Use only commercially available cleaning agents such as mild detergent or alcohol (70% ethanol) to clean the StrepMan. Do not autoclave.
Dimensions	W 13,6 × D 9,7 × H 12 cm; Weight 900 g
Warranty	24 months

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.