



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

Buffer XT-R

Strep-Tactin®XT regeneration buffer

Cat. No.: 2-1045-250

Version: 2.1
Revision Date: 28.07.2021

| | |
|----------------------|---|
| Description | Buffer XT-R is a Strep-Tactin®XT Regeneration Buffer for competitive displacement of biotin from immobilized Strep-Tactin®XT. |
| Components | 3 M MgCl ₂ |
| Form | Solution |
| Concentration | Ready-to-use |
| Amount | 250 ml |
| Stability | 12 months after shipping |
| Storage | 2-8 °C |
| Shipping | Room temperature |
| Hazards | Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided. |
| Application | The buffer is ready to use and should not be diluted. |

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.