

IBA Lifesciences GmbH Rudolf-Wissell-Str. 28 37079 Goettingen Germany Tel.: +49 (0) 551-5 06 72-0 E-mail: info@iba-lifesciences.com www.iba-lifesciences.com

# **Data Sheet**

# pASK-IBA7C

Cat. No.: 2-1326-000

Version: 11.0 Revision Date: 11.06.2021

Description	Expression plasmid. The expression cassette is under transcriptional control of the tetracycline promoter/operator. The expressed recombinant protein will be localized in the cytoplasm.			
Affinity tag	Strep-tag®II fused to the N-terminus of the recombinant protein and can be removed by cleavage with Factor Xa.			
Bacterial Expression	Expression is induced upon addition of 200 $\mu$ g anhydrotetracycline per 1 liter <i>E. coli</i> shaking culture (A <sub>550</sub> = 0.5).			
Expression strain	Any <i>E. coli</i> strain. The <i>tet</i> -promoter works independently from the genetic background of <i>E. coli</i> .			
Resistance	Chloramphenicol			
	<b>Note:</b> The Cam <sup>R</sup> resistance gene codes for homotetrameric chloramphenicol acetyltransferase (MW of the monomer = 26.6 kDa) which is predominantly expressed in the cytosol of <i>E. coli</i> transformed with this plasmid.			
Form	5 μg, dissolved in 20 μl TE buffer, pH 8.0: 10 mM Tris/HCl, 1 mM EDTA			
Concentration	250 ng/μl			
Stability	12 months after shipping			
Storage	recommended: 2-8 °C for frequent usage, -20 °C for long-term storage			
Shipping	room temperature			
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.			

#### For research use only


#### Trademark information

The owners of trademarks marked by "<sup>®</sup>" or "TM" are identified at <u>http://www.iba-lifesciences.com/patents.html</u>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

#### Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <a href="http://www.iba-lifesciences.com/patents.html">http://www.iba-lifesciences.com/patents.html</a> or upon inquiry at <a href="http://info@iba-lifesciences.com/patents.html">info@iba-lifesciences.com/patents.html</a> or upon inquiry at <a href="http://info@iba-lifesciences.com/patents.html">info@iba-lifesciences.com</a> or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.

### Multiple Cloning Site of pASK-IBA7C



**Please note:** Restriction enzymes in bold cut twice. The *Bsa*l sites (isoschizomer of *Eco31*I) at each end of the multiple cloning site are useful for precise and oriented insertion of the recombinant gene by one cleavage reaction only. The "link" contains a restriction site which can be used for subcloning.

### **Features of pASK-IBA7C**

	from bp	to bp
promoter	37	72
forward primer binding site	57	76
Strep-tag <sup>®</sup> II	139	171
Factor Xa cleavage site	172	183
multiple cloning site	184	260
reverse primer binding site	328	344
f1 origin	357	795
CamR resistance gene	917	1576
Tet-repressor	1589	2212
ColE1 origin	2365	2953

Cloning pri	mers for the precise cloning using <i>Bsa</i> l or <i>Eco31</i> I	Sequencing primers:
Forward:	(N <sub>20</sub> ) 5'- NNNNNNGGTCTCNGC GCC NNN NNN	Forward: 5'- GAGTTATTTTACCACTCCCT -3'
Reverse:	$(N_{20})$ 5'- NNNNNNGGTCTCNTA TCA NNN NNN	Reverse: 5'- CGCAGTAGCGGTAAACG -3'