



IBA Lifesciences GmbH
Rudolf-Wissell-Str. 28
37079 Goettingen
Germany
Tel.: +49 (0) 551-5 06 72-0
E-mail: info@iba-lifesciences.com
www.iba-lifesciences.com

Data Sheet

MEXi-TM

Transfection medium for MEXi-293E cells

Cat. No.: 2-6011-010

Version: 3.0
Revision Date: 03.06.2021

Description	MEXi-TM is a transfection medium especially developed for transient transfection and growth of MEXi-293E cells.
Medium Properties	Serum and biotin--free, chemical defined
Endotoxin	<10 EU/ml
Form	Suspension
Amount	1000 ml
Stability	6 months after shipping
Storage	Recommended: 2-8 °C, protect from light
Shipping	Room temperature
Hazards	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
Application	It is applied during the transfection of MEXi-293E cells to obtain the best transfection efficiency.

For research use only

Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.

Important licensing information

This product is covered by intellectual property (IP) rights and on completion of the sale IBA Lifesciences grants respective Limited Use Label Licenses to purchaser. IP rights and Limited Use Label Licenses for said technology are further described and identified at <http://www.iba-lifesciences.com/patents.html> or upon inquiry at info@iba-lifesciences.com or at IBA Lifesciences GmbH, Rudolf-Wissell-Str. 28, 37079 Goettingen, Germany. By use of this product the purchaser accepts the terms and conditions of all applicable Limited Use Label Licenses.