



**IBA Lifesciences GmbH**  
Rudolf-Wissell-Str. 28  
37079 Goettingen  
Germany  
Tel.: +49 (0) 551-5 06 72-0  
E-mail: [info@iba-lifesciences.com](mailto:info@iba-lifesciences.com)  
[www.iba-lifesciences.com](http://www.iba-lifesciences.com)

## Data Sheet

### Buffer CI (10x)

Cat. No.: 6-6320-025, 6-6320-085

Version: 1.4  
Revision Date: 10.05.2024

<b>Description</b>	10x PBS containing 10 mM EDTA and 5% BSA, sterile filtered.
<b>Stability</b>	12 months after shipping; handle under <b>sterile</b> conditions.
<b>Storage</b>	2-8 °C; may also be stored at 15-25 °C
<b>Shipping</b>	Room temperature or blue ice
<b>Hazards</b>	Product is not classified as hazardous according to (EC) No 1272/2008 [CLP]. A Material Safety Data Sheet is provided.
<b>Application</b>	10x Buffer CI is suitable for our Strep-tag® technology approaches for cell isolation and staining.
<b>Usage</b>	Dilute 1:10 with ddH <sub>2</sub> O for 1x Buffer CI working solution.

#### For research use only

##### Trademark information

The owners of trademarks marked by “®” or “TM” are identified at <http://www.iba-lifesciences.com/patents.html>. Registered names, trademarks, etc. used in this document, even when not specifically marked as such, are not to be considered unprotected by law.